REFRIGERATOR

ELETOM79°



KEY FEATURES



User-Interactive & Informative LCD Display



Fast & Precise Cooling Operations



21 CFR Compliance Software with PC Interface



Dual Door Design for Temperature Retention



Email Notifications For Alarms



Efficient Low-Noise Compressor with CFC Free Refrigerant



Process & Event Alarms



High Grade PUF Insulation

FEATURES CONTD...

- ▶ In-built Programmable Timer
- User Programmable Defrost Cycle
- ▶ Intelligent Door-Open Sensing with Alarm
- ▶ Robust Built Quality & Compact Foot-Print
- ≥128 X 64 STN Monochrome Graphic Display
- External Port for Alarm & Acknowledgement
- Very User-Friendly and Informative User Interface
- Comprehensive Alarm System (Process Deviations, Door Open)
- > High Quality PUF Insulation for Better Temperature Retention
- > GMP Model Available



				T	Т	T
Au SERIES	MODELS	EARL-120	EARL-210	EARL-320	EARL-500	EARL-700
DIMENSIONS	CAPACITY (L)	120	210	320	500	700
	INNER CHAMBER (W x D x H) mm	450 x 450 x 600	600 x 600 x 600	600 x 600 x 900	650 x 650 x 1200	750 x 700 x 1350
ACCESSORIES	TRAY	2	2	3	3	4
POWER SUPPLY*	SINGLE PHASE 230 VAC, 50Hz	600 W	750 W	900 W	1200 W	1500 W
TEMPERATURE SENSOR	RANGE	2°C to 8°C				
	RESOLUTION	0.1°C				
	ACCURACY	± 1°C				
CONTROLLER	PID CONTROLLER	Intelligent Microprocessor Based, Tuned PID Temperature Controller				
COOLING	COMPRESSOR	Hermetically Sealed, Low Noise Compressor, with CFC Free Refrigerant used for Efficient Operation				
AIR CIRCULATE	BLOWER	Air Circulated Vertically for Uniform Distribution of Temperature				
ILLUMINATION	LIGHTS	2 LED Lights Available for Illumination of Chamber				
INSTRUMENT MOC	INNER CHAMBER	Inner Chamber Made of SS304 Grade				
	OUTER CHAMBER	Outer Chamber Made of Powder Coated CRCA Sheet				
	OUTER DOOR	Outer Side Made of Powder Coated CRCA Sheet & Inner Side Made of SS304				
	INNER DOOR	Inner Door Made of Superior Quality Toughen Glass for Viewing without Disturbing Environment				
	INSULATION	Insulation Between Inner and Outer Chamber is Filled with High Grade 65mm PUF				
	GASKET	Silicon Based Gasket is Provided for the Door to Avoid Air Leakage				
	VALIDATION PORT	Validation Port is Provided on the Side of the Instrument with Sealed Cap				
	TRAYS	Mesh Type Removable Stainless Steel Trays with User Adjustable Spacing Between Trays				
SAFETY	OPERATIONAL SAFETY	Safety Door Sensor for Operational Cut-Off to Retain Temperature Inside Chamber				
	ELECTRICAL	Electro-Magnetic Interference Filter, Leakage Current Protector, Fast Acting Fuse, MCB				
	SAFETY	Incorporated for User Safety & Component Protection				
	TEMP. SAFETY	Factory Set, Auto Electrical Cut-Off Thermostat Provided				
ACCESSORIES	DOCUMENTS	IQ,OQ Document, NABL Calibration Certificates* for Sensor & Controller				

^{*}For Certain Sensors & Controllers we Provide Non-NABL Certificates.



^{*}Actual Power Values Depends upon Site Voltage Ratings.